

Pollution lab 4

MOHAMED SABBAR



DETERMINATION OF FREE CO₂ IN WATER

Background Information

- ▶ **Carbon dioxide (CO₂)** is present naturally as a result of animal respiration, the decay of organic of matter, and the decomposition certain minerals. It is the major source of acidity in polluted water samples.


Background Information

- ▶ **Surface waters** typically contain less than **10 ppm** (mg/L) dissolved (CO_2), while **ground waters** may contain **several hundred ppm** (mg/L). All aquatic organisms release (CO_2) into the water. Some of it bubbles to the surface, some of it dissolves (mixes in) with the water, but most of the carbon dioxide found in the water is **produced by** organisms (bacteria mostly) that carry on decomposition of dead material.

Background Information

- ▶ During daylight all plants use carbon dioxide and give off oxygen By photosynthesis process. At night, the opposite is Plants use oxygen and give off carbon dioxide. By true respiration process high level of carbon dioxide usually **indicates** that there is a lot of dead material undergoing decomposition. This may occur naturally, but could be the result of different types of water pollution. (CO_2) and water united to form carbonic acid (H_2CO_3), which analyses to bicarbonate ion and positive hydrogen ion as in the formula:

Background Information



Background Information

- ▶ The relationship between pH of the aquatic medium and CO₂ forms can be explained as following:

Medium	PH	CO ₂ Forms
Acidic	Low	Dissolved free form carbonic acid (H ₂ CO ₃)
Neutral	Moderate	Bicarbonate forms (HCO ₃)
Alkaline	High	Carbonate form (CO ₃)

Test Procedure:

- ▶ 1. Collect your water sample carefully to prevent gas volatilization.
- ▶ 2. Take 100 ml from water sample either supplement or irrigated water by cylinder and put it in a flask.
- ▶ 3. Add 10 drops from **phenolphthalein** as indicator solution.
- ▶ 4. Titrate with **0.025N sodium hydroxide** solution (NaOH).
- ▶ 5. Stir the water sample gently during the titration.
- ▶ 6. The point of titration is the start of **pink color** appearance in the solution end
- ▶ 7. Calculate Volume of sample (ml)

$$\text{CO}_2(\text{PPM}) = \frac{\text{ml tirtrant (NaOH)} * (0,025)N * 1000 \text{ mg/L}}{\text{Volume of sample (ml)}}$$