



**Ministry of higher Education &
Scientific Research
Al-Rasheed University College/
Pharmacy Department**



**Practical Inorganic pharmaceutical chemistry I
Third class / 1st semester
(2022-2023)**


**Lab 6
Assay of Borax**

**Done by:
Assist. Lecturer: Jessica sh. Hanna**

Assay of Borax

$\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$ (M.Wt. 381.4)

- Powdered borax is white, consisting of soft crystals that dissolve easily in water. Considered as (weak base)
- Borax, also known as sodium borate, sodium tetraborate, is an important boron compound, and a salt of boric acid.
- Borax synthesized from the reaction of boric acid with sodium carbonate




Pharmaceutical application of borax:

1. Used as buffering agent, and antiseptic.
2. Borax is used as natural preservative
3. For the treatment of fungal foot diseases

Chemical principle:

- Borax samples may sometimes be contaminated with boric acid or sodium carbonate, Thus **two titrations** are carried out;
 1. Borax with Hydrochloric acid (sodium carbonate if present).
 2. Borax with sodium Hydroxide (Boric acid if present).
- Borax is a salt derived from a weak acid and a strong base, so its aqueous solution can be assayed using a standard 0.5N hydrochloric acid solution in (acid-base titration). Titrate Borax with 0.5N HCl Boric acid is liberated



Borax

Boric acid

- The reaction between borax and Hydrochloric acid liberates boric acid a very weak acid that needs to increase its acidity (by converting it into a strong acid) so that it can react with sodium hydroxide by direct titration. This is achieved by the addition of **excess?** of **mannitol** (polyhydroxylated alcohol).

- Mannityl boric acid complex is formed which is strong enough to titrate directly with 1N NaOH.





Figure (1-1) structure of Mannityl boric acid complex

Procedure:

1. dissolve the sample of Borax, in 40 ml. distilled water, add 2 drops of methyl red as indicator.
2. Titrate with 0.5 N HCl
The color changes from **yellow** → **pink**.
3. Boil? and cool the solution, then Record the volume of HCl used.
4. Add 4gm. of mannitol using phenolphthalein as indicator (add 2 drops). and titrate with 1N sodium hydroxide, The color changes from **Pink** → **colorless** → **pink**
Record the volume of NaOH solution used.





Figure (1-1) Titration apparatus.

➤ **Note**

The solution of the titration should be boiled after the endpoint to expel (if any) carbon dioxide generated (if sodium carbonate is present):



This is to prevent the formation of **carbonic acid** which may affect the result of the second titration with sodium hydroxide solution.

Calculation:

Each 1 eq.wt. of NaOH \equiv 1 eq.wt. of Boric acid

Each 4 eq. wt. of boric acid \equiv 1 M. Wt. of Borax

[Each 4 eq. wt. of NaOH \equiv 1 M. Wt. of Borax] /4

Each 1 eq. wt. of NaOH \equiv 1/4 M. Wt. of Borax (**N= no. of eq.wt / V.L**)

V.L * 1N NaOH \equiv 381.4gm * 1/4 of Borax

[1000 ml 1N NaOH \equiv 381.4gm * 1/4 of Borax] /1000

1ml. of 1N NaOH \equiv 381.4

1000*4

Each 1ml. of 1N NaOH equivalent to **0.09535 g.** of Borax.

V_1 = Volume of 0.5N HCl

V_2 = Volume of 1N NaOH

✚ The correct volumes of the two titrations can then tell the presence and type of impurities, besides the quantities of each.

1. If $V_1 = V_2$

This means that the sample of borax is **pure**. Sometimes a difference of **not more than 0.3ml** between the two volumes is allowed, in this case use the average.

Each 1ml. of 1N NaOH equivalent to **0.09535 g. of Borax.**

Wt. of borax = $V \times$ chemical factor

2. If $V_1 > V_2$

The volume of 0.5N HCl solution consumed is more than the volume of 1N NaOH solution consumed. Then the sample of borax is **impure** and contains, in addition to borax, **sodium carbonate**.

$V_1 - V_2$ = the volume of 0.5N HCl solution consumed by sodium carbonate.

Each 1ml. of 1N NaOH equivalent to **0.09535 g. of Borax.**

Each 1ml. of 0.5N HCl equivalent to **0.0265 g. of sodium carbonate.**

3. If $V_2 > V_1$

The volume of 1N NaOH solution consumed is more than the volume of 0.5N HCl solution consumed

Then the sample of borax is **impure** and contains in addition to borax, **boric acid**

$V_2 - V_1$ = the volume of 1N NaOH solution consumed by the contaminant boric acid.

Each 1ml. of 0.5N HCl equivalent to **0.09535 g. of Borax**

Each 1ml. of 1N NaOH equivalent to **0.06184 g. of boric acid.**