Experiment No.(5) Extraction

- **Extraction** is a method of separation and purification of organic compounds that depends on the ability of the compound to dissolve in two immiscible solvents, e.g. water and chloroform.

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
Partition coefficient: *is the ratio of the organic solvent concentration in the two layers at a certain temperature.*

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A good solvent for extraction should:

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
1. Separatory funnel.
2. Erlenmeyer flasks.

Procedure:

1. Ensure that the stopcock at the bottom of the separatory funnel is closed.
2. Add solvent (5ml of CCl_4 which contains a low percentage of acetic acid as an impurity) to the separatory funnel.
3. Add (15ml) of sodium bicarbonate solution with concentration of (10%) to the separatory funnel.
4. Insert the stopper to top of the separatory funnel.
5. Remove the funnel from the ring and grasp its neck with one hand, holding the stopper down firmly. Invert funnel and

immediately open stopcock to release pressure, close the stopcock, and mix the layers by shaking the funnel.

6. Place the separatory funnel in the ring once more and wait until the layers have completely separated. Remove the stopper and open the stopcock to draw off the bottom layer (which is the organic layer containing CCl_4) into a labeled Erlenmeyer flask.
7. Steps to the experimental procedure are shown in (figure 1).


Figure 1 Extraction experimental procedure.