

EXPERIMENT 6
OXIDATION-REDUCTION REACTION
(KMNO₄ STANDARDIZATION)

Introduction

- These types of titration depend upon the change in oxidation number of reacting substances, a redox reaction between an oxidizing agent and a reduction agent.
- A substance that is oxidized is termed a *reducing agent*, while the substance that is reduced is term an *oxidizing agent*.

Oxidation

Is defined as a loss of one or more electrons by atom or an increase in the oxidation number.

Reduction

Is defined as a gain of one or more electrons by atom or a decrease in the oxidation number (opposite oxidation).

Oxidation can never take place without reduction, because something must be able to pick up the electrons lost by the oxidized molecules.

Types of Sterilizer

- Many *sterilizer* have the property of killing bacteria because they are either oxidizing or reducing agent.

Oxidizing agent

Chlorine

Potassium permanganate

Iodine solution

Reducing agent

Formaldehyde

Sulfur dioxide

Procedure

- KMnO_4 is a strong oxidizing agent, and will be standardized in this experiment versus oxalic acid using sodium oxalate as the primary standard.

- **Procedure:**

1. Take 0.13 gm of sodium oxalate to 250 ml conical flask.
 2. Add 20 ml of 1 N sulfuric acid (H_2SO_4).
 3. Heat the solution to about 60 °C.
 4. Then titration solution versus permanganate slowly.
- The solution must still be hot at the end point.



Calculation

■ Calculate the normality of KMnO_4

$$\text{■ Titer factor of } \text{NaOX} / \text{KMnO}_4 = \frac{\text{Wt. NaOX}}{\text{Vol. KMnO}_4}$$

$$\text{■ Normality of } \text{KMnO}_4 = \frac{T \times 1000}{\text{GEW}}$$