



Al-Rasheed University College
Pharmacy Department
2nd Stage / 2nd Course
(2019 – 2020)



TOTAL WBC COUNTING

PHYSIOLOGY LAB.
2ND STAGE / 2ND COURSE
LAB. (4)

INTRODUCTION:

- Leukocytes are originated in bone marrow from stem cells (hemocytoblasts)
- Once they are fully mature, they leave bone marrow and enter the blood stream
- Normal value: 4000 – 11000 WBCs / Imm^3 of blood
- When talking about WBCs, there are important terms that must be known:
 - Leukocytosis
 - Leukopenia
 - Leukemia

INTRODUCTION (CONTINUE):

■ **Leukocytosis:**

- It means an increase in the number of WBCs above upper normal limit
- Mostly occurs in
 1. acute infections (severe bacterial infection, severe malarial infection)
 2. After hemorrhage
 3. During pregnancy
 4. Postoperative (after surgery)
 5. And some types of anemia
- Total WBCs count up to 20000 WBCs / 1mm^3 of blood is considered leukocytosis

INTRODUCTION (CONTINUE):

■ **Leukopenia:**

- A decrease in the number WBCs below the lower normal limit
- Mostly occurs in
 1. Viral infection
 2. Typhoid fever
 3. After radiation therapy and chemotherapy
 4. Chronic use of steroid medication
 5. Pernicious anemia
- Total WBCs count below 4000 WBCs / Imm^3 of blood is considered leukopenia

INTRODUCTION (CONTINUE):

■ **Leukemia:**

- Is characterized by uncontrolled increase and proliferation of one type or more of the different hematopoietic cells, and associated with many changes in the circulating cells of the blood
- The most 2 common types of leukemia are:
 - Lymphocytic leukemia
 - Myelocytic leukemia
- Leukemia can be further subdivided into:
 - Acute: disease progressive rapidly and there is marked morphological changes in cells
 - Chronic: the disease is slow and there is no obvious changes
- Total WBCs count exceeding 20000 WBCs / Imm^3 of blood could be considered leukemia

IMPORTANT NOTES:



- Total WBCs counting can be done manually using (special slide, stain and microscope)
- and can be done by automated method (hematology analyzer)
- Manual method in total WBCs counting is still superior in certain circumstances like (counting WBCs in cerebrospinal fluid and synovial fluid)

OBJECTIVE (AIM) OF THE EXPERIMENT:

1. The presence of infection and follow the progress of certain diseases
2. Presence of Leukemia
3. Monitor the body response to various treatments
4. Conditions or diseases that may weaken the immune system, like HIV infection, chemotherapy, radiation therapy causing a drop in WBC count.
5. Monitor bone marrow function


BASIC REQUIREMENTS:

- 70% alcohol and cotton
- Sterile blood lancet
- Special WBC pipette
- Special slide (hemocytometer)
- Cover slip
- Special dilution colored fluid (Turk's solution)
- Microscope



PROCEDURE:

I- WBC PIPETTE



- Small and graduated (0.5 IU – 1.1 IU)
- Used to withdraw blood sample, and then dilute it with **Turk's solution**

PROCEDURE:

2-TURK'S SOLUTION



- A special dilution colored fluid
- Composed of:
 - Glacial acetic acid (for destruction of RBCs)
 - Gentian violet or methylene blue (for staining WBCs)
 - Distilled water (for dilution of sample)



PROCEDURE:

3-HEMOCYTOMETER (COUNTING CHAMBER)

- Hemocytometer is a special slide containing a counting chamber
- It is used for manual measurement of different types of cells in blood and other body fluids
- Used to determine the total number of:
 - WBCs in a blood sample, cerebrospinal fluid sample, and synovial sample
 - RBCs in a blood sample
 - Platelets in a blood sample
 - Sperms in a seminal fluid sample



HOW TO CALCULATE THE TOTAL NUMBER OF WBC?

- **Part one: derived from the (WBC pipette)**

0.5 unit of blood is diluted to 11 units with Turk's solution

- 0.5 unit (B) \rightarrow 11 unit (B+T)

Mix by rotating the pipette for 2 – 5 min then **discards 3 drops which are equivalent to 1 (IU) unit.**

- 0.5 unit (B) \rightarrow 11 unit (B+T) – 1 unit (B+T)
- 0.5 unit (B) \rightarrow 10 unit (B+T)
- $[0.5 \text{ unit (B)} \equiv 10 \text{ unit (B+T)}] \times 2$
- 1 unit (B) \equiv 20 unit (B+T)

Since $[1 \text{ mm}^3 \text{ (B)} \equiv 20 \text{ mm}^3 \text{ (B+T)}] / 20$

So $1/20 \text{ mm}^3 \text{ (B)} \equiv 1 \text{ mm}^3 \text{ (B+T)}$

$$1 \text{ mm}^3 \text{ (B+T)} \equiv 1/20 \text{ mm}^3 \text{ (B)}$$



HOW TO CALCULATE THE TOTAL NUMBER OF WBC?

■ Part Two: derived from (hemocytometer)

Calculate the volume of (1) corner square

$$\begin{aligned} V &= L(\text{length}) * W(\text{width}) * D(\text{depth}) \\ &= 1\text{mm} * 1\text{mm} * 0.1\text{mm} \\ &= 0.1 \text{ mm}^3 (B+T) \end{aligned}$$

<u>Volume</u>	<u>Cells</u>
0.1 mm ³ (B+T)	W
1 mm ³ (B+T)	?

Then 1 mm³ (B+T) = 10W


Since 1 mm³ (B+T) \equiv 1/20 mm³ (B) ----- (from the pipette equation)



Then:


1/20 mm ³ (B)	10W
1 mm ³ (B)	?

$$W1+W2+W3+W4 / 4 = \text{Result (Average W)}$$


$$1 \text{ mm}^3 (B) = 200W$$







— count
--- don't count



EXAMPLE:

- Calculate the total WBC count for a blood sample when the total count for square A= 80, square B= 100, square C= 120 and square D= 88 cells. State if the result is normal, leukocytosis, leukopenia or leukemia?

Answer:

$$A + B + C + D = ?$$

$$80 + 100 + 120 + 88 = 388$$

$$388 / 4 = 97 \text{ (average)}$$

$$\text{Since } 1\text{mm}^3 \text{ (B)} = 200 \text{ W}$$

$$\text{Then } 200 * 97 = 19400 \text{ WBC} / 1\text{mm}^3 \text{ of blood}$$

(Leukocytosis) because the normal value is

$$4000 - 11000 \text{ WBC} / 1\text{mm}^3 \text{ of blood}$$

EXAMPLE:

- calculate the total WBC if the total WBCs of two chambers (A and C) is 48 and state if the result is normal, leukocytosis, leukopenia or leukemia?

Answer: we have 2 squares

$$A + C = 48$$

$$48 / 2 = 24 \text{ (average)}$$

$$\text{Since } 1 \text{ mm}^3 \text{ (B)} = 200 \text{ W}$$

$$\text{Then } 200 * 24 = 4800 \text{ WBC} / 1 \text{ mm}^3 \text{ of blood}$$

the result is **normal** because the normal value is

$$4000 - 11000 \text{ WBC} / 1 \text{ mm}^3 \text{ of blood}$$

EXAMPLE:

- calculate the total WBC count for a blood sample when the total count of square A= 19, square B= 20, and square C= 10 cells. State if the result is normal, leukocytosis, leukopenia or leukemia?

Answer 1: we have 3 squares

We will choose (A and C) squares only

$$A + C = ?$$

$$19 + 10 = 29$$

$$29 / 2 = 14.5 \text{ (average)}$$

Since $1\text{mm}^3 \text{ (B)} = 200 \text{ W}$

Then $200 * 14.5 = 2900 \text{ WBC} / 1\text{mm}^3 \text{ of blood}$

Result: *Leukopenia* because normal value is
4000 – 11000 WBC/ 1mm^3 of blood

Answer 2: we have 3 squares (the best answer)

We will take the summation of all (3) squares

$$A + B + C = ?$$

$$19 + 20 + 10 = 49$$

$$49 / 3 = 16.3 \text{ (average)}$$

Since $1\text{mm}^3 \text{ (B)} = 200 \text{ W}$

Then $200 * 16.3 = 3260 \text{ WBC} / 1\text{mm}^3 \text{ of blood}$

Result: *leukopenia* because normal value is
4000 – 11000 WBC/ 1mm^3 of blood



STAY SAFE!

Keep Your Distance!

Thank you
for your
patience

